

INTRODUCING

The Newest Innovation for OR Efficiency



Folds down to easily fit the sterilization tray



The innovative design of the Osher Sterile Field Magnifier allows the technician to inspect fine surgical instrument tips prior to handing off to the surgeon during the procedure. With its large 50 mm lens, the magnifier is also ideal for loading IOLs into the injector and aligning the irrigation sleeve on the phaco or I/A tip prior to use. Manufactured from stainless steel and fine quality optical glass, the magnifier can be cleaned and sterilized by conventional means and easily fits into most surgical instrument trays for ease of processing.



Use the Osher Sterile Field Magnifier to:

- Inspect microscopic instrument tips for consistency and damage prior to use
- Assist OR staff in differentiating microsurgical instrument tips
- Easily and accurately load IOLs into the injector system with its powerful lens
- Accurately align and adjust the irrigation sleeve on the phaco tip or on the IA tip
- View diamond blades for cleanliness and edge consistency
- Can also be used at the instrument cleaning station for inspection during processing
- Base dimensions: 65 mm x 140 mm (2.6" x 5.5")

Order your Osher Sterile Field Magnifier now!

Item #

E7002

Description

Osher Sterile Field Magnifier